

Personal Protection

Apparel

Bouffant Caps

Caps are made of spunbonded polypropylene. They are lightweight, cool, and comfortable with a soft elastic headband.

Surgeons Caps

Caps are constructed of a non-woven and absorbent material with a spunbonded polypropylene top for comfort and breathability. Traditional tie on style. Fluid resistant.

Face Masks

Molded-style **Cone Masks** are fluid-resistant and have an adjustable aluminum nose piece.

Earloop Procedure Masks feature a pleated style.

Earloop Procedure Masks with Anti-fog Shield provide a clear plastic wrap-around that protects eyes from fluid.

Tie-On Surgical Masks features a pleated style with two ties top and bottom.

Isolation Gowns

Comfortable, fluid-resistant protection. Full length coverage, elastic arm cuffs, and ties at neck and waist

Impervious Gowns

Polyethylene gown with thumb loop is impervious to fluids.

Shoe Covers

Spunbonded polypropylene shoe covers with sewn seams. Elastic opening provides a snug fit. Non-skid soles offer improved floor traction.





Personal Protection

Apparel

Pro Advantage has you covered from head-to-toe. The personal protection apparel products feature quality and performance. These single-use products are designed for comfort and prevention of cross-contamination. The entire product line is latex-free.

Item #	Description	Size	Color	UOM
BOUFFANT CAPS				
P701010	Bouffant Cap	21"	White	100/bx, 5 bx/cs
P701015	Bouffant Cap	21"	Blue	100/bx, 5 bx/cs
P701020	Bouffant Cap	24"	White	100/bx, 5 bx/cs
P701025	Bouffant Cap	24"	Blue	100/bx, 5 bx/cs
SURGEON	IS CAPS			
P701030	Surgeons Cap		Blue	100/bx, 10 bx/cs
FACE MA	SKS			
P703010	Molded-Style Cone Mask		Blue	50/bx, 6 bx/cs
P703020	Earloop Procedure Mask		Blue	50/bx, 10 bx/cs
P703025	Earloop Procedure N			, , , , ,
	with Anti-Fog Shield		Blue	25/bx, 4 bx/cs
P703015	Tie-On Surgical Mask		Blue	50/bx, 12 bx/cs
ISOLATIO	N GOWNS			
P704030	Isolation Gown	Regular	Blue	10/bg, 5 bg/cs
P704020	Isolation Gown	Regular	Yellow	10/bg, 5 bg/cs
P704025	Isolation Gown	X-Large	Yellow	10/bg, 5 bg/cs
IMPERVIOUS GOWNS				
P704110	Impervious Gown	Regular	Blue	15/bg, 5 bg/cs
P704115	Impervious Gown	X-Large	Blue	15/bg, 5 bg/cs
SHOE CO	VERS			
P702010	Shoe Covers	Universal	Blue	100/bx, 3 bx/cs
P702015	Non-Skid	Universal	Blue	100/bx, 3 bx/cs
P702020	Non-Skid	X-Large	Blue	100/bx, 3 bx/cs
Distributed by:				















